

A high-angle, close-up photograph of three surgeons in an operating room. They are wearing blue surgical gowns, masks, and caps. The scene is dimly lit with a strong blue light source, creating a clinical and focused atmosphere. The surgeons are huddled together, looking down at a patient who is not visible in the frame.

IRRAS

Cranial Access Kit

IRRAflow[®]

• *Hummingbird[®]*
ICP Monitoring

A kit designed to better meet your cranial access needs.

IRRAS CRANIAL ACCESS KIT IC40-D (WITH DRUGS)

Preparation Components

- 2 Disposable Razor
- 1 DuraPrep Solution, 6ml*
- 10 Gauze Sponges, 4" x 4", 12-ply
- 2 Medicine Cup, 2oz, 60cc
- 2 4" x 4" 6-ply Gauze
- 1 Providone Iodine, 1oz Foil Pouch
- 2 Providone Iodine Swabs
- 2 Culture Tube, 16x125mm

Cranial Access Preparation

- 1 Ruler 6"
- 1 Marking Pen
- 1 Fenestrated Drape, 24" x 26", Clear
- 1 Half Drape, 53" x 77"
- 2 CSR Wrap, 36" x 36"
- 3 Absorbent Towel, 17" x 20"

Cranial Access

- 1 Lidocaine with EPI, 30ml Ampule*
- 2 Sodium Chloride, 10ml Ampule*
- 2 Syringe, 12cc
- 2 Needle, 18 G x 1 1/2"
- 1 Needle, 25 G x 5/8"
- 1 Spinal Needle, 18 G x 3 1/2"
- 1 Scalpel, #15
- 1 Scalpel, #11
- 1 Retractor, Blunt
- 1 Hand Crank Drill
- 1 4mm Drill Bit with Stop and Hex Wrench
- 1 7.1mm Drill Bit with Stop and Hex Wrench
- 1 Mosquito Forceps, Curved

Wound Closing/Dressing

- 1 Bone Wax
- 1 Tegaderm
- 1 2.0 Silk Suture*
- 1 3.0 Nylon Suture*
- 1 Needle Holder, Serrated
- 1 Adson Forceps
- 1 Adson Forceps with Teeth
- 1 Suture Scissors

ORDERING INFORMATION

Product #	Product Name	Description
IC40-D	IRRAS Cranial Access Kit IC40-D (with Drugs)	Cranial Access Kit with Sodium Chloride and Lidocaine
IC40	IRRAS Cranial Access Kit (without Drugs)	Cranial Access Kit without Sodium Chloride and Lidocaine

IRRAS

www.iras.com



Customer Service
1-800-946-0458
US.customerservice@iras.com

IRRAS
11975 El Camino Real, Suite 304
San Diego, CA 92130
Tel: (800) 213-4604

*Specific components included in kit IC40-D only
HM-101_A